



System One

Category: CPAP

REMstar Auto



milnadumedicalsystem.com

Optimizing CPAP prescription pressures at home

System One features

- Advanced event detection*
- Opti-Start feature
- A-Flex for comfort
- System One Resistance Control
- Mask fit check
- Waveform reports
- Oximetry module
- On-board memory
- Smart ramp
- Encore patient management system
- SD card and wired, wireless, and broadband modem capable

Humidification

- System One humidification control
- Pre-heat humidifier option
- Heated Tube humidifier option

REMstar Auto advantages

- CPAP, Auto, CPAP-Check, Auto-Trial modes
- A-Flex for comfort
- Opti-Start technology
- Advanced event detection
- Modem and SD card capable
- Optional, integrated oximetry available
- Heated Tube option

*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).

Home Healthcare
Respiratory Solutions

★ Company of the Year ★
2012-2014

awarded by
FROST & SULLIVAN

PHILIPS


RESPIRONICS

REMstar Auto specifications

| Part number | Description |
|-------------|---|
| N561S | REMstar Auto A-Flex with SD card |
| IN561HS | REMstar Auto A-Flex with humidifier and SC card |
| IN561TS | REMstar Auto A-Flex with Heated Tube humidifier and SD card |

| Specifications | |
|---------------------------------|--|
| Pressure range | 4-20 cm H ₂ O |
| Flex pressure relief | 0 to 3 |
| Ramp time | 0 to 45 min (5-min increments) |
| Starting ramp pressure | 4 to CPAP; patient adjustable (CPAP mode) 4 to min pressure (Auto mode) |
| Dimensions | 18 cm L x 14 cm W x 10 cm H |
| Weight | 1.36 kg (without humidifier) |
| Filters | Foam and optional ultra-fine |
| Device set-up | LCD/control wheel/push button |
| Data storage capacity (minimum) | Display: compliance 30-day averages SD card > 1 year On-board: up to 6 months |
| Compliance meter | Breathing detection |
| Altitude compensation | Automatic |
| Electrical requirements | 100-240 VAC, 50/60 Hz |
| Humidification | Integrates with System One 60 Series heated humidifier and System One 60 Series Heated Tube humidifier |
| Warranty | 2 years |
| Miscellaneous | Lighted LEDs |

| Sleep therapy accessories | |
|---|-------------|
| Humidification | Part number |
| System One 60 Series humidifier | IN6H |
| System One 60 Series Heated Tube humidifier | IN6T |
| Water chamber (includes both 50 and 60 Series lids) | 1063785 |
| Heated Tube humidifier upgrade kit | 6T-CONVKIT |
| Tubing | Part number |
| System One Heated Tube | SYS1HT15 |
| System One Performance Tube (15mm) | 1074400 |
| Performance Tubing (22mm) | 1032907 |
| Filters | Part number |
| Pollen (foam) filter, reusable (1/pack) | 1035443 |
| Ultra-fine filter, disposable (1/pack) | 1063096 |
| Accessories | Part number |
| Device side panel | 1099592 |
| Carrying case | 1063857 |
| Power cord Europe | 1039014 |
| Power cord UK | 1039026 |



Respironics, A-Flex, Auto-Trial, C-Flex, CPAP-Check, Encore, REMstar, and System One are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/respironics



© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Broudy NB 12/19/13 MCI 4105940 PN 1113946

Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics Benelux
+33 1 47 28 30 83
Philips Respironics Deutschland
+49 8152 93 06 0
Philips Respironics France
+33 2 51 89 36 00
Philips Respironics Iberia
+34 91 566 90 00

Philips Respironics Italy
+39 039 203 1
Philips Respironics Sweden
+46 8 120 45 900
Philips Respironics Switzerland
+41 6 27 45 17 50
Philips Respironics United Kingdom
+44 800 1300 845
Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668 USA