



System One

Category: CPAP

REMstar Pro



milnadumedicalsystem.com

Your long-term CPAP solution

System One features

- Advanced event detection*
- C-Flex+ for comfort
- System One Resistance Control
- Mask fit check
- Waveform reports
- Oximetry module
- On-board memory
- Encore patient management system
- SD card and wired, wireless, and broadband modem capable

Humidification

- System One humidification control
- Pre-heat humidifier option
- Heated Tube humidifier option

REMstar Pro advantages

- CPAP, CPAP-Check
- C-Flex+ for comfort
- Advanced event detection
- Modem and SD card capable
- Optional, integrated oximetry available
- Heated Tube option

*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).



As your Allies in Better Sleep and Breathing, we lead the way by providing more than advanced technology: we provide exceptional support.

PHILIPS

RESPIRONICS

REMstar Pro specifications

Part number	Description
IN461S	REMstar Pro C-Flex+ with SD card
IN461HS	REMstar Pro C-Flex+ with humidifier and SD card
IN461TS	REMstar Pro C-Flex+ with Heated Tube humidifier and SD card

Specifications	
Pressure range	4-20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	4 to CPAP; patient adjustable (CPAP mode) 4 to min pressure (Auto mode)
Dimensions	18 cm L x 14 cm W x 10 cm H
Weight	1.36 kg (without humidifier)
Filters	Foam and optional ultrafine
Device set-up	LCD/control wheel/push button
Data storage capacity (minimum)	Display: compliance 30-day averages SD card > 1 year On-board: up to 6 months
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100-240 VAC, 50/60 Hz
Humidification	Integrates with System One 60 Series heated humidifier and System One 60 Series Heated Tube humidifier
Warranty	2 years
Miscellaneous	Lighted LEDs

Sleep therapy accessories	
Humidification	Part number
System One 60 Series humidifier	IN6H
System One 60 Series Heated Tube humidifier	IN6T
Water chamber (includes both 50 and 60 Series lids)	1063785
Heated Tube humidifier upgrade kit	6T-CONVKIT
Tubing	Part number
System One Heated Tube	SYS1HT15
System One Performance Tube (15mm)	1074400
Performance Tubing (22mm)	1032907
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1035443
Ultrafine filter, disposable (1/pack)	1063096
Accessories	Part number
Device side panel	1099592
Carrying case	1063857
Power cord Europe	1039014
Power cord UK	1039026

Respironics, A-Flex, Auto-Trial, C-Flex, CPAP-Check, Encore, REMstar, and System One are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/respironics



© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Broudy NB 4/1/14 MCI 4105939 PN 1113945

Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics Benelux
+33 1 47 28 30 83
Philips Respironics Deutschland
+49 8152 93 06 0
Philips Respironics France
+33 2 51 89 36 00
Philips Respironics Iberia
+34 91 566 90 00

Philips Respironics Italy
+39 039 203 1
Philips Respironics Sweden
+46 8 120 45 900
Philips Respironics Switzerland
+41 6 27 45 17 50
Philips Respironics United Kingdom
+44 800 1300 845
Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668 USA