

Home Medix[®]

Oxygen Concentrator



- Ultra Quiet Technology (36dbA)
- Nebulizer Function
- Purity Sensor
- Intelligent Diagnosis System
- Time Setting Control

Models : HM - 5

HM - 8

HM - 10



**Best Choice for
Home Oxygen**

Home Medix[®] Oxygen Concentrator

Specifications

Model	HM - 5	HM - 8	HM - 10
Flow Rate	1 - 5 L/M	1 - 8 L/M	1 - 10 L/M
Concentration	93 ± 3%		
Power	390 W	450 W	610 W
Power Supply	220 ± 10% - 50Hz ± 1%		
Sound Level	≤ 36 db	≤ 50 db	≤ 50 db
Outlet Pressure	45 kpa ± 10%		
Net Weight	21.5 kg	24 kg	24 kg
Gross Weight	23.5 kg	26 kg	26 kg
Size	372 x 340 x 612 mm		
Warranty	3 Years		
Colour	Creamy		
Alarm	Power failure Alarm; Temperature Alarm; Low Purity Alarm		
Optional Functions	SPO ₂ Sensor		
Accessories	1x Power Line, 1x Inlet Filter Cotton, 1x High Efficiency Filter, 1x Nebulizer, 2x Nasal Pipe, 1x Oxygen Connecting Tube, 1x Humidification Bottle, 1x User manual, 1x Warranty Card.		

*Specifications subject to change without notice



Marketed By :

TAMILNADU MEDICAL SYSTEM



Head Office : No.11/26, G3, Swathy Flats, CRM Street, Varadarajapuram, Ambattur, Chennai - 600 053.

Branch Office : No. 6/31, Kadappa Road, New Lakshmiapuram, Kolathur, Chennai - 600 099.

Customer Care No : 93809 31624 / 93805 26938

E-Mail : taminadumedicalsystem@gmail.com

Web : www.tamilnadumedicalsystem.com